



TOR's / Specifications of Dialysis Kits

01	<p>Hemodialyzer</p> <p>From macro-design to membrane: Nanostructure; the case of the FX-class of hemodialyzers</p> <p>Detail: Hollow Fiber Dialyzer Polysulfone (Steam Sterilized) with Tubing (0.8mm,1.0mm,1.3mm,1.4mm,1.6mm,1.8mm)with 6.5mm/6.3 tubing</p>
02	<p>Heparin Sodium Inj</p> <p>Indication:</p> <ul style="list-style-type: none"> • Prophylaxis of thrombo-embolism • Use as anticoagulant in the therapy of acute venous arterial thromboembolism • Prevention of blood clotting during use of extracorporeal circulation (heart-lung machine, haemodialysis) <p>Features:</p> <ul style="list-style-type: none"> • Contains benzyl alcohol as preservative • Aqueous solution <p>Presentation: 5 ml vial Inj. Heparin 5ml</p> <p>Product Description:</p> <ul style="list-style-type: none"> • Product Form: Solution • Pack Size: 4ltr • Marketed By: MEDIPAK PHARMA • Generic Category: Dialysis Solution • Ingredients: Acid & Basic Concentrate for Bicarbonate Haemo Dialysis <p>Bicarbonate Dialysis Solution part A+B</p>
03	<p>Biocarbonate</p> <p>Features:</p> <ul style="list-style-type: none"> • Medical grade PVC Tubing. • Various diameter and type for different machines is available. • Pump Segment for Conventional and high flow rates. • It is easy to handle and kink-resistant tubing. • Non-DEHP, Non-pyrogenic and Non-toxic. <p>Specification:</p> <ul style="list-style-type: none"> • Pump Segment: 6.6x9.8mm & 8x12mm • Chamber: 22 / 30mm • Main Tube: 4.6x6.8mm • Accessory: Pillow, Transducer Protector, IV Set, Drainage Bag – one way, two way • Standard Packing: 24sets /carton <p>Make: Imported</p>
04	<p>Blood Line</p> <p>Features</p> <p>Ultra-thin needle walls for painless insertion and superior blood flow</p>
05	<p>AV Needle Standard</p>


 Dr. Faruk Shariq
 Chief Medical Officer (CFO) North
 Oil & Gas Development Co. Ltd
 Islamabad

		<p>Textured wings provide secure grip during needle manipulation</p> <p>Soft flexible tubing absorbs shock and reduces needle movement in patient</p> <p>Color coded wings for easy recognition of gauge size</p> <p>Easy to use occlusion clamp</p> <p>Luer tapered for easier connection with all bloodlines</p> <p>AV Fistula Needle with Clamps size 16G,17G,18G OR equivalent</p>
--	--	--


 Dr. Medicinal Officer (F/P) North
 Oil & Gas Development Co. Ltd
 Oil City, Islamabad